

Designed by a Gynecologist for Gynecologists!

VALTCHEV®

UTERINE MOBILIZER

MODEL VUM - 6



VUM-6 Sterilization Tray

ST-6



Product Features:

- Interior is clearly marked with catalogue numbers, simplifies counting & re-ordering
- Made of top quality Radel® plastic
- Able to be sterilized in an autoclave
- Easy to use

Product Features:

- Allows 135° of anteversion and 15° of retroversion
- Ergonomic design allows for single-handed use and improved line-of-sight
- Reusable, cost-effective solution eliminates need for disposable instruments
- Includes a compatible Valtchev® Tenaculum Forcep
- Includes four uterine obturator sizes and two sizes of cannulas

Contact Conkin Surgical for Valtchev® parts, repairs and regular servicing.

Email: sales@conkinsurgical.com

Phone: 604.638.2757